



JOLI Diagnostic Inc · 1348 North Forest Road · Williamsville · NY
14221 · (716) 639-0443

tTG Human Tissue Transglutaminase IgG

Method:

IgG class antibodies can be detected in patient serum using standard EIA methods.

Specimen Requirements and Shipping:

1. Approximately 0.5 mL of patient serum is necessary to perform the analysis. The serum should be placed in a small, tightly-capped tube. The tube should be labeled with the patient name, date of birth, ID number, date of specimen collection and referring physician name.
2. The sample should be kept frozen and stored at -20°C or -70°C until it can be shipped on dry ice to our facility.
3. Package according to regulations regarding shipment of infectious material and pack enough dry ice to last 48 hours.
4. Ship overnight by appropriate courier (i.e. UPS/Fed Ex) to the address above.

Please note: The laboratory is closed on weekends and holidays. When shipping, please ship overnight and ensure a **weekday** delivery.

Turn-Around Time: Test Performed Weekly

Samples are batched and analyzed weekly. The results will be faxed and/or mailed to you immediately following assay completion.

Billing is made to the institution from which the sample was sent.
The cost of this analysis is \$45.00.